

1
2
3
4
5
6
7
8
9
10
11
12
13
14
15
16
17
18
19
20
21
22
23

**Target Product Profile for Misoprostol and Mifepristone Point-of-Care
Quality Testing**

Disease Area: Medical Abortion

Intervention: (near) point-of-care quality detection devices

Version: 1.2 March 2026

24 **Table of Contents**

25 VERSION: 1.2 MARCH 2026.....1

26 **1 BACKGROUND3**

27 1.1 General3

28 1.2 Purpose of this Target Product Profile4

29 **2 SUMMARY: INTERVENTION INDICATION AND TARGET USERS5**

30 **3. TPP 1 (UPPER LEVEL OF SUPPLY CHAIN): CORE VARIABLES7**

31 **4. TPP 2 (LOWER LEVEL OF SUPPLY CHAIN): CORE VARIABLES26**

32 **3 REFERENCES43**

33

DRAFT

34 1 Background

35 1.1 General

36 Comprehensive abortion care relies on access to quality-assured, affordable medical abortion medicines —
37 mifepristone, misoprostol, and co-packaged mifepristone-misoprostol (combi-pack) — including for self-managed
38 abortion. Increasing the availability of these medicines is essential to improve access and depends on medicines
39 policies and regulation, procurement and distribution, as well as health worker and end-user knowledge. Ensuring
40 the quality of abortion medicines, such as mifepristone and misoprostol, is essential for their safe and efficacious
41 use.

42 Both quality assurance and quality control play crucial roles in this process, though they serve distinct functions.
43 Quality assurance is a proactive, process-driven approach that ensures abortion medicines are manufactured,
44 stored, and distributed under strict regulatory standards. It includes compliance with Good Manufacturing
45 Practices (GMP), supplier qualification, and continuous monitoring of production systems to prevent quality issues
46 before they arise. Quality assurance helps ensure that only safe and efficacious products enter the supply chain.
47 Quality control, in contrast, is a product-focused process that involves testing medicines at different points in the
48 supply chain to verify they meet established standards. This includes analyzing samples for active pharmaceutical
49 ingredient (API) content, detecting impurities, and assessing stability. Quality control is critical to confirming that
50 mifepristone and misoprostol tablets meet quality standards at the time they reach the end-user.

51 Abortion medicines are distributed through diverse channels, including pharmacies, clinics, telemedicine services
52 (via post and internet services), and community health programs. In some settings, they are legally regulated,
53 while in others, access is restricted due to legal and social barriers. To ensure that good quality medical abortion
54 medicines are used, there is a need to advance quality testing methods and technologies. The current approach
55 to conducting quality surveys is to sample medicines in the field and send those samples to a laboratory for quality
56 testing. Several studies have been conducted in the last decade, identifying substandard misoprostol and
57 mifepristone in multiple countries, which can lead to inadequate efficacy (lower amount of API leading to reduced
58 effectiveness) and safety concerns (such as impurities leading to increased risk of adverse reactions)¹⁻³.
59 Stigmatized products such as medical abortion medicines are also more likely to be falsified. Indeed, falsified
60 products that do not contain any API have been found in circulation and remain a critical concern⁴. However, all
61 these quality survey studies are time-consuming and costly, using high-performance liquid chromatography, and
62 do not readily inform local decision-making, especially at the point of care. In addition, in some regulatory settings,
63 resources may be limited and a national or certified laboratory may not be available to test suspected substandard
64 medicines. While samples can sometimes be sent to laboratories in neighboring countries, obtaining results may
65 take considerable time; in such situations, a quantitative device capable of providing numerical values of the API
66 content in a tablet could offer valuable preliminary information. Having tools for quality control at the point of
67 care would be valuable for real-time monitoring that enables prompt action and decision-making to remove
68 substandard and falsified products from the supply chain. However, there are currently no methods for testing
69 mifepristone and misoprostol at this level.

70 Point-of-care (POC) or near POC technologies have been developed for a range of medicines, including infectious
71 disease medicines, and could be applied to abortion medicines^{5,6}. Currently, fingerprint technologies are available
72 to identify whether a given API is present in a finished pharmaceutical product (FPP) and spectroscopy can be
73 used to evaluate an FPP by comparing it to a reference formulation. These technologies, however, are qualitative
74 and do not quantify the API content of the product, a critical parameter to ensure product quality, adequate
75 dosing and, ultimately, efficacy.

76

77 1.2 Purpose of this Target Product Profile

78 Target product profiles (TPPs) specify the minimum and optimal requirements for new health products, including
79 diagnostic tests⁷⁻¹¹. They are used by WHO and other international organizations to address major public health
80 gaps, including the need to spur development and accelerate implementation of health products⁷⁻¹¹. TPPs are an
81 important resource for funders, researchers, product developers, manufacturers and regulators¹². They guide key
82 stakeholders on the characteristics required of new health products to meet pre-specified clinical and public
83 health needs. They inform research and development strategies, help frame product dossiers, streamline
84 communication with regulatory agencies and help funders set targets¹³. In addition, medical products approved
85 by the US Food and Drug Administration that addressed pre-specified TPPs have been linked to more rapid
86 regulatory review¹⁴.

87
88 There are currently no TPPs available for POC or near POC technologies that test the quality of medical abortion
89 medicines. Availability of POC devices that can test the quality of medical abortion tablets would facilitate the
90 implementation of safer and more effective access to abortion care. Development of this TPP provides clear and
91 consistent requirements to help drive innovation, research, and implementation of an accessible quality testing
92 device for medical abortion medicines that meet global health needs.

93
94

DRAFT

95 2 Summary: Intervention Indication and Target Users

96

97 This TPP focuses on quality testing of **misoprostol** as a priority use-case because it is more widely used and
98 distributed than mifepristone, including for medical abortion and postpartum hemorrhage management.
99 Misoprostol is also chemically less stable, degrading quickly when exposed to heat and humidity, and its low active
100 ingredient content (measured in micrograms) makes accurate quality testing technically challenging. In contrast,
101 mifepristone is more stable and typically manufactured and distributed through more controlled channels.
102 Therefore, developing a robust, field-appropriate quality testing method for misoprostol addresses a greater
103 public health need and it would pave the way for mifepristone quality testing in the near future.

104

105 In developing a TPP for a medical abortion product testing device, we recognize that users at different points along
106 the supply chain have different needs and product requirements. We have therefore created two separate TPPs.

107

108 **TPP 1 (Upper level of supply chain)** – for use by individuals involved in the procurement and wholesale distribution
109 of medical abortion products (e.g. Ministries of Health, county governments, large distributors, and other
110 centralized supply agencies). In addition, the device can be used by regulatory authorities for field screening.

111

112 *Upper level of supply chain scenario:*

113 When a shipment of misoprostol tablets is received, stakeholders at the upper level of the supply chain can sample
114 and test tablets on-site with this device, or it can be used at storage facilities to measure the quality before further
115 distribution. If the device indicates that a batch contains substandard levels of the active ingredient, the
116 stakeholders can take immediate action: the suspect batch can be quarantined or rejected outright instead of
117 being distributed. The stakeholder may then invoke quality safeguards – for example, renegotiating with the
118 manufacturer or supplier based on the test findings, or enforcing contract clauses that require acceptable quality
119 before payment is finalized. Such proactive measures strengthen accountability in the supply chain by signaling to
120 manufacturers that substandard products will not be tolerated. Furthermore, any batch flagged as substandard
121 by the device should be withheld from downstream distribution.

122 While this device is not a replacement for compendial laboratory testing, its findings can serve as an early warning:
123 regulatory bodies can be notified to investigate and perform formal analytical testing on the suspected products.

124

125

126 **TPP 2 (Lower level of supply chain)** – for use by individuals involved at point of distribution or
127 use of medical abortion products (e.g. pharmacies, hospitals, community centers).

128

129 *Lower level of supply chain scenario:*

130 When a pharmacy, clinic, or community health center receives misoprostol tablets—especially from informal or
131 non-centralized sources—the provider can use this device to test a sample tablet on-site before dispensing it to
132 patients. If the device indicates that the tablet contains insufficient active ingredient, the provider can
133 immediately set aside the suspect stock and avoid selling or using it. They may also contact the distributor to
134 report the issue, request a replacement or refund, and choose to source from a different supplier in the future. In
135 some cases, the result may prompt the provider to stop working with a supplier altogether. This enables frontline
136 health workers to play an active role in ensuring medicine quality, even in settings where regulatory oversight is
137 weak or formal quality assurance is lacking.

138

139 The device can also be used as a spot-check tool to verify the quality of products already in stock, or to build trust
140 with patients concerned about counterfeit or ineffective medicines. While the tool is not a substitute for official
141 regulatory testing, it empowers providers to detect substandard products early and take action to protect patient
142 safety.

143
144 Upper-level supply chain settings (ministries of health, large distributors, etc.) generally have more budget, more
145 controlled environments, better technical infrastructure and trained staff. By contrast, lower-level supply chain
146 settings (pharmacies, clinics or community health care centers) need a low-cost, portable, easy-to-use screening
147 tool with minimal training. Regulatory and distribution contexts might also differ: in some countries medical
148 abortion medicines are strictly channeled through official systems, while in others they circulate informally.
149 Because quality problems are common, having two TPPs lets us tailor requirements for each level – such as
150 balancing cost, training level, and portability for centralized procurement sites/ labs versus point-of-care settings.

151
152 The testing device is intended for use by various stakeholders across the supply chain to assess the quality of
153 medical abortion medicines. It does not replace the formal analytical testing required for regulatory approval or
154 official quality certification, but supports flagging potential substandard or falsified products.

DRAFT

155 **3. TPP 1 (Upper level of supply chain): Core Variables**

156

Variable	Minimum	Optimum	Annotations
Product Overview & Target Context			
<p>1: Intended use</p>	<p>A device to assess the presence of sufficient active pharmaceutical ingredient (API) of misoprostol in a tablet at the near point-of-care.</p>	<p>A device to assess the presence of sufficient active pharmaceutical ingredient of misoprostol and mifepristone in a tablet at the near point-of-care.</p>	<p>In practical use, when a shipment of misoprostol tablets is received, stakeholders at the upper level of the supply chain can sample and test tablets on-site with this device.</p> <p>If the device indicates that a batch contains substandard levels of the active ingredient, the stakeholders can take immediate action: the suspect batch can be quarantined or rejected outright instead of being distributed. The stakeholder may then invoke quality safeguards – for example, renegotiating with the manufacturer or supplier based on the test findings, or enforcing contract clauses that require acceptable quality before payment is finalized.</p> <p>Such proactive measures strengthen accountability in the supply chain by signaling</p>

Variable	Minimum	Optimum	Annotations
			<p>to manufacturers that substandard products will not be tolerated.</p> <p>Furthermore, any batch flagged as substandard by the device should be withheld from downstream distribution.</p> <p>While this device is not a replacement for compendial laboratory testing, its findings can serve as an early warning: regulatory bodies can be notified to investigate and perform formal analytical testing on the suspected products.</p> <p>Letrozole is excluded from the list of target products prioritized for measurement in this TPP because additional evidence is needed on its safety, effectiveness, and acceptability.</p>
<p>2: Target Use Setting</p>	<p>Upper level of the supply chain, including procurement and wholesale distribution settings of medical abortion products (e.g. Ministries of Health, county governments, large distributors, and other centralized supply agencies). In addition, the device can be used by regulatory authorities for field screening.</p>	<p>Same as minimum</p>	<p>The tablets can be tested before being distributed to the abortion care providers (healthcare facilities, pharmacies)</p>

Variable	Minimum	Optimum	Annotations
<p>3:</p> <p>Target operator and skills</p>	<p>Individuals working at upper-level supply chain settings, including the procurement level/ wholesale distributor’s location or associated laboratories</p> <p>Requires a basic skill set and training depending on the background of person and technology</p>	<p>Same as minimum</p> <p>Requires minimal skill set and training, depending on background of person and technology</p>	<p>The responsibility of quality testing cannot be shifted onto users themselves, especially when dealing with sensitive medications such as abortion drugs, where individuals may not have the means, knowledge, or safe conditions to verify quality. Users must be able to trust that the medicines they obtain are genuine and effective, without facing the added burden of testing or uncertainty.</p> <ul style="list-style-type: none"> • No skills: Lay person with no healthcare or technical background • Minimal skills: Lay person or staff member after brief demonstration • Basic skills: Trained pharmacy or healthcare support staff. • Advanced skills: Healthcare professional • Specialist skills: Laboratory-qualified professional
<p>Product Design & Functionality</p>			

Variable	Minimum	Optimum	Annotations
<p>4: Portability</p>	<p>Moveable by one person, a small table-top device that can be carried/moved with or without the help of a cart.</p>	<p>Handheld device.</p>	<p>Handheld is the optimistic option for maximum access.</p> <p>Minimum target for the benchtop device is based on the current weight and dimension of a mobile mini-laboratory¹⁵, and airline allowed weight for check-in luggage. Optimistic weight has been based on the maximum recommended weight for a woman to manually lift¹⁶.</p> <p>Minimum target for the handheld device is based on the current weight and dimension of a Handheld Raman Analyzer¹⁷. Optimistic weight and dimensions are based on a smartphone.</p>
<p>5: Required accessories</p>	<p>None</p>	<p>None</p>	<p>Reducing dependence on companion parts simplifies deployment, especially in rural areas with limited supply chains or import restrictions.</p> <p>Additional accessories could include for example a probe, sample holder, or UV viewing box. These accessories do not relate to those needed for sample preparation</p>
<p>6: Training requirements</p>	<p>Basic operation (instrument handling, sample preparation, running measurements): maximum 2 day of hands-</p>	<p>Basic operation (instrument handling, sample preparation, running measurements): maximum 1 day of remote or self-training (for example</p>	

Variable	Minimum	Optimum	Annotations
	<p>on training is sufficient for someone with prior lab experience.</p> <p>Data interpretation and routine quality control use: requires maximum 2 weeks to get comfortable with</p>	<p>online video instructions) is sufficient for someone with prior lab experience.</p> <p>Data interpretation and routine quality control use: requires maximum 1 week to get comfortable with</p>	
<p>7: Languages required for operating</p>	<p>English Language of target countries</p>	<p>Same as minimum</p>	
<p>8: Results display and interpretation</p>	<p>Requires companion device provided by the user (no onboard display or result output)</p> <p>Result is displayed immediately e.g. as yes/no or as % API detected. Results are easy and intuitive to interpret.</p>	<p>Integrated display with onboard result interpretation or companion device is not needed for result display.</p> <p>Result is displayed immediately e.g. as yes/no or as % API detected. Results are easy and intuitive to interpret.</p>	
<p>9: Connectivity</p>	<p>Test results can be transferred to a laptop/ tablet using a cable or by using Bluetooth.</p>	<p>Same is minimum</p>	<p>There is a need to transfer results to device (computer/ tablet/smartphone), in order to ensure adequate storage of the test results and printing if needed.</p> <p>The output should be universal, e.g. USB3</p>

Variable	Minimum	Optimum	Annotations
10: Internal results storage	100 tests	500 tests	Requirement on the minimum numbers of results that can be stored within the device before it runs out of storage capacity. *Only applicable for multi use devices
11: Data security	Reference libraries, algorithms and generated data should be appropriately protected.	Same as minimum	*Only applicable for multi use devices
12: Software - Interoperability standards and format	Software runs on commonly available operating systems including their older versions. Data output is given in several appropriate file formats, including .csv. Software updates are provided free of charge	Device can be operated by a smartphone or no software for analysis of results required. Data output is given in several appropriate file formats, including .csv. Software updates are provided free of charge and are available for older versions of commonly available operating systems.	*Only applicable if software is required.
13: Software functionality & user access	Single-user application with basic interface; fixed data output. No differentiation between user types. No automatic data sharing through a centralized dashboard.	Multi-user platform with role-based access Includes data sharing between (inter-) national users, centralized dashboard to view aggregated results, and automatic notifications/alerts when substandard products are detected.	*Only applicable if software is required. For example, users at higher-level settings can access detailed analytical data, and exportable reports (e.g., PDFs for data filing), while lower-level users receive simplified outputs (e.g., pass/fail or semi-quantitative).

Variable	Minimum	Optimum	Annotations
14: Support	Full support through a website chat box, WhatsApp/Wechat assistance, with email replies guaranteed within 24 hours	Full support through a website chat box, WhatsApp/Wechat assistance, and a toll-free number, with email replies guaranteed within 24 hours	
15: Sample preparation	<p>Minimal destruction of tablet sample is needed for the device to measure the quality.</p> <p>Tablet must be free from packaging (box and blister)</p> <p>Could require simple sample preparation (e.g. no usage of solvent use besides distilled water or alcohol, maximum three steps).</p> <p>Sample preparation can be done using basic laboratory equipment, such as a mortar, pestle and vortex mixer</p>	<p>No destruction of tablet sample is needed for the device to measure the quality.</p> <p>Tablet must be free from packaging (box only, measurement through blister)</p> <p>No additional sample preparation required.</p>	<p>Misoprostol tablets are typically sealed individually in aluminum blister packs, which protect them from moisture and humidity, as the tablets are highly sensitive to moisture.</p> <p>Primary packaging consists of aluminum blisters, while secondary packaging typically includes a small cardboard box.</p>
16: Reagents Requirement	Reagents might be required, but are easily accessible in LMIC	No reagents are required.	As mentioned under ‘sample preparation’, reagents may include distilled water or alcohol, which in general are available at basic pharmacies
17: Measurement Time	Under 10 minutes	Under 5 minutes	The measurement time does not include sample preparation

Variable	Minimum	Optimum	Annotations
Technology & Performance			
Single use device			
18A: Target product	Device can only be used for misoprostol tablets	Same as minimum Plus: Device can be used for mifepristone tablets	This TPP focuses on misoprostol for the quality testing device as a priority because it is more widely used and distributed than mifepristone, including medical abortion and postpartum hemorrhage management. Misoprostol is also chemically less stable, degrading quickly when exposed to heat and humidity, and its low active ingredient content (measured in micrograms) makes accurate quality testing technically challenging. A single use device is defined as a disposable one time use only
19A: Quantitation	Semi-quantitative	Same as minimum	The device is not intended to replace the formal analytical testing required for regulatory approval or official quality certification. Definitions

Variable	Minimum	Optimum	Annotations
			<ul style="list-style-type: none"> • Semi-quantitative: provides a yes/no or range rather than an exact amount. • Quantitative: provides an exact numerical measurement of the API amount present. <p>For a single-use device, operation and interpretation should be as simple as possible, giving a clear binary result:</p> <ul style="list-style-type: none"> • Yes — tablet contains sufficient misoprostol API. • No — tablet contains insufficient misoprostol API.
20A: Analytical specificity/ selectivity	No interference from other excipients commonly used in the target product (misoprostol tablets or mifepristone tablets) formulations.	Same as minimum	
21A: Limit of Detection (LOD)	90% of the label claim	90% of the label claim	Misoprostol tablets for medical abortion are labeled to contain 200 micrograms of API. International quality standards require 90–110% of the label claim, corresponding to 180–220 micrograms. → Tablets containing less than 180 micrograms (90%) of misoprostol API are considered substandard products.

Variable	Minimum	Optimum	Annotations
			<p>The Limit of Detection (LOD) is the lowest amount that a device can detect.</p> <p>For example, considering the label claim of misoprostol is 200 micrograms per tablet:</p> <ul style="list-style-type: none"> • 90% LOD: 180 micrograms is the lowest amount that the device can detect. • 50% LOD: 100 micrograms is the lowest amount that the device can detect. • 10% LOD: 20 micrograms is the lowest amount that the device can detect. <p>To enable a semi-quantitative device that indicates whether a misoprostol tablet is of good quality (using the cut-off value of 180 micrograms), the LOD should be set at 90%.</p>
<p>22A: Limit of Quantification (LOQ)</p>	<p>Not applicable for semi-quantitative devices</p>	<p>Not applicable for semi-quantitative devices</p>	<p>*Only applicable for quantitative tests</p>
<p>23A: Sensitivity and specificity for detecting the target</p>	<p>Sensitivity 80 % Specificity 80 %</p>	<p>Sensitivity 95% Specificity 95%</p>	<p>A minimum sensitivity and specificity of ≥80% is acceptable because it reflects a reasonable balance between technical feasibility and public health benefit—allowing the device to identify most severely substandard tablets</p>

Variable	Minimum	Optimum	Annotations
substance at the limit of detection			<p>(well below the 90–110% label claim) while recognizing that some margin of error is tolerable in low-resource settings. This level is consistent with established benchmarks for triage or first-line tools in global health, where perfect accuracy is not always achievable but early detection still provides significant value.</p> <p>At the optimal level, ≥95% sensitivity and specificity ensures the device can reliably distinguish between very low-API tablets and those within or near the acceptable range. This higher threshold supports confident decision-making, minimizes false positives and negatives, and aligns with standards used for high-quality diagnostics intended to trigger regulatory or supply chain action.</p>
24A: Calibration	Calibration of the device is not needed	Calibration of the device is not needed	
25A: Robustness	Method should be robust to different analytical procedure parameters, ambient temperature, humidity and duration of the procedure.	Same as minimum	
26A: Reproducibility	Method should give reproducible results when device is used by a different person,	Same as minimum	

Variable	Minimum	Optimum	Annotations
	on a different day and in a different place after transport.		
27A: Maintenance Requirements	No maintenance needed	Same as minimum	
28A: Internal Control/ Reference	Optional internal control/ reference items (e.g. values, standards) are either included or inexpensive and readily available in LMIC.	No internal control or reference items needed	Built-in controls improve reliability, ensure test validity, and enhance user trust— particularly critical when devices are used in decentralized settings.
29A: Companion equipment	No companion equipment is required	No companion equipment is required	Eliminating reliance on external equipment simplifies implementation and improves usability in clinics with limited infrastructure and in home settings.
Multi use device			
18B: Target product	Device can only be used for misoprostol tablets Device has the potential to measure other APIs	Same as minimum Plus: Device can be used for mifepristone tablets Device can be used for testing other API, for which methodologies and reference libraries exist (or can be developed)	For higher end technologies, such as the NIR (Near-Infrared) or Raman technologies, methods could be developed and validated to measure other API’s using the same device.

Variable	Minimum	Optimum	Annotations
19B: Quantitation	Semi-quantitative	Quantitative	<p>The device is not intended to replace the formal analytical testing required for regulatory approval or official quality certification.</p> <p>A quantitative version of the device could be useful at upper levels settings, where users may include Ministries of Health, procurement agencies, and regulators. In some settings, access to accredited laboratory facilities may be limited or unavailable, meaning tablets cannot easily be sent for confirmatory testing. In such situations, a more precise numerical result can support evidence-based decisions on product acceptance, distribution, investigation of suspected substandard medicines, and supply chain quality monitoring.</p> <p>Definitions</p> <ul style="list-style-type: none"> • Semi-quantitative: provides a yes/no or range rather than an exact amount. • Quantitative: provides an exact numerical measurement of the API amount present.
20B:	No interference from other excipients commonly used in the target product	Same as minimum	

Variable	Minimum	Optimum	Annotations
Analytical specificity/ selectivity	(misoprostol tablets or mifepristone tablets) formulations.		
21B Limit of Detection (LOD)	90% of the label claim	10% of the label claim (assuming a quantitative device)	<p>Misoprostol tablets for medical abortion are labeled to contain 200 micrograms of API. International quality standards require 90–110% of the label claim, corresponding to 180–220 micrograms. → Tablets containing less than 180 micrograms (90%) of misoprostol API are considered substandard products.</p> <p>The Limit of Detection (LOD) is the lowest amount that a device can detect (provide yes/no). The Limit of Quantification (LOQ) is the lowest amount that a device can measure (provide a number)</p> <p>For example, considering the label claim of misoprostol is 200 micrograms per tablet:</p> <ul style="list-style-type: none"> • 90% LOD/LOQ: 180 micrograms is the lowest amount that the device can detect (LOD) or measure (LOQ). • 50% LOD/LOQ: 100 micrograms is the lowest amount that the device can detect (LOD) or measure (LOQ).

Variable	Minimum	Optimum	Annotations
			<ul style="list-style-type: none"> 10% LOD/LOQ: 20 micrograms is the lowest amount that the device can detect (LOD) or measure (LOQ). <p>To enable a semi-quantitative device that indicates whether a misoprostol tablet is of good quality (using a cut-off value of 180 micrograms), the LOD should be set at 90%.</p> <p>To support the use case of quantitative device in the upper level settings, the LOQ could be set at 100 micrograms (50% of the label claim). This allows the device to generate reliable numerical results across a clinically and programmatically relevant range, enabling higher-level users to distinguish between compliant products, moderately degraded tablets, and clearly substandard medicines.</p> <p>In general, the LOD is about one tenth of the LOQ. For example, a 50% LOQ (100 micrograms) corresponds to a 5% LOD (10 micrograms). However, because misoprostol is present at very low microgram levels, a 10% LOQ is considered more feasible.</p>
22B:	Not applicable for semi-quantitative devices	50% of the label claim	<p>*Only applicable for quantitative tests</p> <p>See explanation at 'Limit of Detection (LOD)' variable</p>

Variable	Minimum	Optimum	Annotations
Limit of Quantification (LOQ)			
23B: Sensitivity and specificity for detecting the target substance at the limit of detection	Sensitivity 80% Specificity 80%	Sensitivity 95% Specificity 95%	<p>A minimum sensitivity and specificity of $\geq 80\%$ is acceptable because it reflects a reasonable balance between technical feasibility and public health benefit—allowing the device to identify most severely substandard tablets (well below the 90–110% label claim) while recognizing that some margin of error is tolerable in low-resource settings. This level is consistent with established benchmarks for triage or first-line tools in global health, where perfect accuracy is not always achievable but early detection still provides significant value.</p> <p>At the optimal level, $\geq 95\%$ sensitivity and specificity ensures the device can reliably distinguish between very low-API tablets and those within or near the acceptable range. This higher threshold supports confident decision-making, minimizes false positives and negatives, and aligns with standards used for high-quality diagnostics intended to trigger regulatory or supply chain action.</p>
24B: Calibration	Calibration of the device is needed, can be performed by the user	Calibration of the device is needed, can be performed by the user	

Variable	Minimum	Optimum	Annotations
	Calibration time is maximum 10 minutes	Calibration time is maximum 5 minutes	
25B: Robustness	Method should be robust to different analytical procedure parameters, ambient temperature, humidity and duration of the procedure.	Same as minimum	
26B: Reproducibility	Method should give reproducible results when device is used by a different person, on a different day and in a different place after transport.	Same as minimum	
27B: Maintenance Requirements	LMIC-based third party vendors can provide maintenance. Spare parts are available within 4 weeks after placement of the order Recalibration or servicing are needed every two years	Maintenance can be performed by user itself Spare parts are available within 2 weeks after placement of the order Recalibration or servicing are needed every three years	Devices requiring minimal maintenance but with accessible support pathways are essential for long-term use in resource-limited settings. Reliance on international service may not be sustainable or responsive.
28B: Internal Control/ Reference	Optional internal control/ reference items (e.g. reference library, values, standards) are either included or inexpensive and readily available in LMIC.	No internal control needed	Built-in controls improve reliability, ensure test validity, and enhance user trust— particularly critical when devices are used in decentralized settings.
29B: Companion equipment	Small, portable, reusable tabletop devices may be required to process samples for tests.	Small, portable, reusable handheld device, or no companion equipment required.	Eliminating reliance on external equipment simplifies implementation and improves usability in clinics with limited infrastructure and in home settings.
Cost & Financial Sustainability			

Variable	Minimum	Optimum	Annotations
31: Device cost	Single use device: Under \$10 Multi use device: Under \$3000	Single use device: Under \$5 Multi use device: Under \$250	Upfront cost of the device/technology. Should be priced appropriately for intended markets, or supplier should be willing to negotiate on pricing.
32: Cost per test	If consumables and/or single use items (for example sample holder, cartridge) are required, one test should cost less than \$10.	If consumables and/or single use items (for example sample holder, cartridge) are required, one test should cost less than \$5. If no consumables and/or single use items are required, there are no additional costs per test	*Excluding costs for personnel
33: Cost of software	Software is offered at an affordable cost with discount options for low-resource settings.	Software is offered free of charge and made available open-source.	*Only applicable if software is required
34: Accessibility	Company provides shipment to every country with appropriate documentation for import of device and consumables.	Same as minimum	
35: Availability	Robust supply chain management for devices and consumables with minimized risk of stockouts in case of increased demand.	Same as minimum	

Variable	Minimum	Optimum	Annotations
	<p>Components to operate the device are readily available from several LMIC sources.</p>		
<p>36: Volume estimates</p>	<p>Volumes compatible with targeted use by central procurement entities in high-need regions.</p> <p>This assumes initial adoption by national Ministries of Health, large distributors, and central supply agencies to test misoprostol quality for medical abortion in key markets, covering the most critical supply channels without extending to every outlet.</p>	<p>Volumes compatible with broad integration across supply chains for misoprostol (used for several indications including medical abortion, PPH, and gastric ulcer management).</p> <p>This scenario reflects uptake by procurement agencies, umbrella health facilities, NGOs, and large pharmacy networks, supported by inclusion in national quality control programs.</p>	<p>For higher-end analytical platforms (e.g., near-infrared spectroscopy), broader interest and purchasing willingness are expected, as these devices could later be extended to test other APIs beyond misoprostol—such as mifepristone, oxytocin, or other essential medicines—once methods are developed. While this may not directly increase initial testing volumes, it enhances value justification and procurement rationale for investing in a versatile, multi-analyte testing system.</p>
<p>37: Stability / Shelf Life / Life cycle</p>	<p>Single use device: Shelf life ≥ 1 years at 30–45°C and 75–90% RH.</p> <p>Multi use device: The lifespan of the device is approximately 5 years, depending on usage intensity, environmental conditions, and maintenance</p>	<p>Single use device: Shelf life ≥ 2 years at 30–45°C and 75–90% RH.</p> <p>Multi use device: The lifespan of the device is approximately 10 years, depending on usage intensity, environmental conditions, and maintenance</p>	<p>For supply chain reliability, products should remain functional across broad temperature/humidity conditions typical of LMICs. Shelf stability without cold chain simplifies logistics and minimizes stockouts.</p> <p>Computers and tablets, potential companion devices, usually start showing issues after 4-5 years, QC instruments have a longer life cycle but need to be replaced after a certain period of use</p>

157

158

159

160 **4. TPP 2 (Lower level of supply chain): Core Variables**

161

162

Variable	Minimum	Optimum	Annotations
Product Overview & Target Context			
<p>1: Intended use</p>	<p>A device to confirm the presence of sufficient active pharmaceutical ingredient (API) of misoprostol in a tablet at the point-of-care.</p>	<p>A device to confirm the presence of sufficient active pharmaceutical ingredient of misoprostol and mifepristone in a tablet at the point-of-care.</p>	<p>In practical use, when a pharmacy, clinic, or community health center receives misoprostol tablets—especially from informal or non-centralized sources—the provider can use this device to test a sample tablet on-site before dispensing it to patients.</p> <p>If the device indicates that the tablet contains insufficient active ingredient, the provider can immediately set aside the suspect stock and avoid selling or using it. They may also contact the distributor to report the issue, request a replacement or refund, and choose to source from a different supplier in the future. In some cases, the result may prompt the provider to stop working with a supplier altogether. This enables frontline health workers to play an active role in ensuring medicine quality,</p>

Variable	Minimum	Optimum	Annotations
			<p>even in settings where regulatory oversight is weak or formal quality assurance is lacking.</p> <p>The device can also be used as a spot-check tool to verify the quality of products already in stock, or to build trust with patients concerned about counterfeit or ineffective medicines.</p> <p>While the tool is not a substitute for official regulatory testing, it empowers providers to detect substandard products early and take action to protect patient safety.</p> <p>Letrozole is excluded from the list of target products prioritized for measurement in this TPP because additional evidence is needed on its safety, effectiveness, and acceptability.</p>
<p>2: Target Use Setting</p>	<p>Setting where the medical abortion medicines are distributed to the buyer, or the place where the medicines are being provided (e.g. pharmacies, hospitals, community centers, mobile points).</p>	<p>Same as minimum</p>	<p>The tablets can be tested before being provided to the buyer either once received at the facility or in real time.</p>
<p>3: Target operator and skills</p>	<p>Individuals working at pharmacies, hospitals, community centers</p>	<p>Same as minimum</p>	<p>The responsibility of quality testing cannot be shifted onto users themselves, especially when dealing with sensitive medications such as abortion drugs, where individuals may not</p>

Variable	Minimum	Optimum	Annotations
	Requires minimal skill set and training, depending on background of person and technology	User with no laboratory, medical, pharmaceutical or other related training can operate the device/test with minimal training using clear instructions (written, visual, or digital)	<p>have the means, knowledge, or safe conditions to verify quality. Users must be able to trust that the medicines they obtain are genuine and effective, without facing the added burden of testing or uncertainty.</p> <ul style="list-style-type: none"> • No skills: Lay person with no healthcare or technical background • Minimal skills: Lay person or staff member after brief demonstration • Basic skills: Trained pharmacy or healthcare support staff. • Advanced skills: Healthcare professional • Specialist skills: Laboratory-qualified professional
Product Design & Functionality			
4: Portability	Handheld device, can be carried in a backpack	Handheld device, can be carried in a backpack	<p>Handheld is the optimistic option for maximum access.</p> <p>Minimum target for the benchtop device is based on the current weight and dimension of a mobile mini-laboratory¹⁵, and airline</p>

Variable	Minimum	Optimum	Annotations
			<p>allowed weight for check-in luggage. Optimistic weight has been based on the maximum recommended weight for a woman to manually lift ¹⁶.</p> <p>Minimum target for the handheld device is based on the current weight and dimension of a Handheld Raman Analyzer ¹⁷. Optimistic weight and dimensions are based on a smartphone.</p>
<p>5: Required accessories</p>	<p>None</p>	<p>None</p>	<p>Reducing dependence on companion parts simplifies deployment, especially in rural areas with limited supply chains or import restrictions.</p> <p>Additional accessories could include for example a probe or sample holder. These accessories do not relate to those needed for sample preparation</p>
<p>6: Training requirements</p>	<p>Basic operation (instrument handling, sample preparation, running measurements): maximum two days of remote or self-training (for example online video instructions) is sufficient for someone without prior lab experience.</p>	<p>Basic operation (instrument handling, sample preparation, running measurements): maximum ½ day of remote or self-training (for example online video instructions) is sufficient for someone without prior lab experience.</p>	

Variable	Minimum	Optimum	Annotations
	Data interpretation and routine quality control use: requires maximum ½ week to get comfortable with	Data interpretation and routine quality control use: requires maximum 1 day to get comfortable with	
7: Languages required for operating	English Language of target countries	Same as minimum	
8: Results display and interpretation	Integrated display with onboard result interpretation or companion device is not needed for result display. Result is displayed immediately e.g. as yes/no. Results are easy and intuitive to interpret.	Same as minimum	
9: Connectivity	Test results can be transferred to a laptop/ tablet using a cable or by using Bluetooth.	Same is minimum	There is a need to transfer results to device (computer/ tablet/smartphone), in order to ensure adequate storage of the test results and printing if needed. The output should be universal, e.g. USB3
10: Internal results storage	50 test	100 tests	Requirement on the minimum numbers of results that can be stored within the device before it runs out of storage capacity. *Only applicable for multi use devices

Variable	Minimum	Optimum	Annotations
11: Data security	Generated data should be appropriately protected.	Same as minimum	*Only applicable for multi use devices
12: Software - Interoperability standards and format	Software runs on commonly available operating systems including their older versions. Data output is given in several appropriate file formats, including .csv. Software updates are provided free of charge	Device can be operated by a smartphone or no software for analysis of results required. Data output is given in several appropriate file formats, including .csv. Software updates are provided free of charge and are available for older versions of commonly available operating systems.	*Only applicable if software is required.
13: Software functionality & user access	Single-user application with basic interface; fixed data output. No differentiation between user types. No automatic data sharing through a centralized dashboard.	Multi-user platform with role-based access Includes data sharing between (inter-) national users, centralized dashboard to view aggregated results, and automatic notifications/alerts when substandard products are detected.	*Only applicable if software is required.
14: Support	Full support through a website chat box, WhatsApp/Wechat assistance, with email replies guaranteed within 24 hours	Full support through a website chat box, WhatsApp/Wechat assistance, and a toll-free number, with email replies guaranteed within 24 hours	

Variable	Minimum	Optimum	Annotations
15: Sample preparation	<p>Minimal destruction of tablet sample is needed for the device to measure the quality.</p> <p>Tablet must be free from packaging (box and blister)</p> <p>No additional sample preparation required.</p>	<p>No destruction of tablet sample, is needed for the device to measure the quality</p> <p>Tablet must be free from packaging (box only, measurement through blister)</p> <p>No additional sample preparation required.</p>	<p>Misoprostol tablets are typically sealed individually in aluminum blister packs, which protect them from moisture and humidity, as the tablets are highly sensitive to moisture.</p> <p>Primary packaging consists of aluminum blisters, while secondary packaging typically includes a small cardboard box</p>
16: Reagents Requirement	No reagents are required.	No reagents are required.	As mentioned under ‘sample preparation’, reagents may include distilled water or alcohol, which in general are available at basic pharmacies
17: Measurement Time	Under 5 minutes	Under 1 minute	The measurement time does not include sample preparation
Technology & Performance			
Single use device			

Variable	Minimum	Optimum	Annotations
18A: Target product	Device can only be used for misoprostol tablets	Same as minimum Plus: Device can be used for mifepristone tablets	This TPP focuses on misoprostol for the quality testing device as a priority because it is more widely used and distributed than mifepristone, including medical abortion and postpartum hemorrhage management. Misoprostol is also chemically less stable, degrading quickly when exposed to heat and humidity, and its low active ingredient content (measured in micrograms) makes accurate quality testing technically challenging. A single use device is defined as a disposable one time use only
19A: Quantitation	Semi-quantitative	Same as minimum	The device is not intended to replace the formal analytical testing required for regulatory approval or official quality certification. Definitions <ul style="list-style-type: none"> • Semi-quantitative: provides a yes/no or range rather than an exact amount. • Quantitative: provides an exact numerical measurement of the API amount present. For a single-use device, operation and

Variable	Minimum	Optimum	Annotations
			interpretation should be as simple as possible, giving a clear binary result: <ul style="list-style-type: none"> • Yes — tablet contains sufficient misoprostol API. • No — tablet contains insufficient misoprostol API.
20A: Analytical specificity/ selectivity	No interference from other excipients commonly used in the target product (misoprostol tablets or mifepristone tablets) formulations.	Same as minimum	
21A: Limit of Detection (LOD)	90% of the label claim	90% of the label claim	Misoprostol tablets for medical abortion are labeled to contain 200 micrograms of API. International quality standards require 90–110% of the label claim, corresponding to 180–220 micrograms. → Tablets containing less than 180 micrograms (90%) of misoprostol API are considered substandard products. The Limit of Detection (LOD) is the lowest amount that a device can detect. For example, considering the label claim of misoprostol is 200 micrograms per tablet:

Variable	Minimum	Optimum	Annotations
			<ul style="list-style-type: none"> • 90% LOD: 180 micrograms is the lowest amount that the device can detect. • 50% LOD: 100 micrograms is the lowest amount that the device can detect. • 10% LOD: 20 micrograms is the lowest amount that the device can detect. <p>To enable a semi-quantitative device that indicates whether a misoprostol tablet is of good quality (using the cut-off value of 180 micrograms), the LOD should be set at 90%.</p>
22A: Limit of Quantification (LOQ)	Not applicable for semi-quantitative devices	Not applicable for semi-quantitative devices	*Only applicable for quantitative tests
23A: Sensitivity and specificity for detecting the target substance at the limit of detection	Sensitivity 80% Specificity 80%	Sensitivity 95% Specificity 95%	

Variable	Minimum	Optimum	Annotations
24A: Calibration	Calibration of the device is not needed	Calibration of the device is not needed	
25A: Robustness	Method should be robust to different analytical procedure parameters, ambient temperature, humidity and duration of the procedure.	Same as minimum	
26A: Reproducibility	Method should give reproducible results when device is used by a different person, on a different day and in a different place after transport.	Same as minimum	
27A: Maintenance Requirements	No maintenance needed	Same as minimum	
28A: Internal Control/ Reference	Optional internal control/ reference items (e.g. values, standards) are either included or inexpensive and readily available in LMIC.	No internal control or reference items needed	Built-in controls improve reliability, ensure test validity, and enhance user trust— particularly critical when devices are used in decentralized settings.
29A: Companion equipment	No companion equipment is required	Same as minimum	Eliminating reliance on external equipment simplifies implementation and improves usability in clinics with limited infrastructure and in home settings.
Multi use device			
18B:	Device can only be used for misoprostol tablets	Same as minimum	For higher end technologies, such as the NIR (Near-Infrared) or Raman technologies,

Variable	Minimum	Optimum	Annotations
Target product	Device has the potential to measure other APIs	Plus: Device can be used for mifepristone tablets Device can be used for testing other API, for which methodologies and reference libraries exist (or can be developed)	methods could be developed and validated to measure other API's using the same device.
19B: Quantitation	Semi-quantitative	Semi-quantitative	The device is not intended to replace the formal analytical testing required for regulatory approval or official quality certification. Definitions <ul style="list-style-type: none"> • Semi-quantitative: provides a yes/no or range rather than an exact amount. • Quantitative: provides an exact numerical measurement of the API amount present. For a lower level setting, assuming users have no or limited training, operation and interpretation should be as simple as possible, giving a clear binary result: <ul style="list-style-type: none"> • Yes — tablet contains sufficient misoprostol API. • No — tablet contains insufficient misoprostol API.

Variable	Minimum	Optimum	Annotations
20B: Analytical specificity/ selectivity	No interference from other excipients commonly used in the target product (misoprostol tablets or mifepristone tablets) formulations.	Same as minimum	
21B: Limit of Detection (LOD)	90% of the label claim	90% of the label claim	The Limit of Detection (LOD) is the lowest amount that a device can detect. For example, considering the label claim of misoprostol is 200 micrograms per tablet: <ul style="list-style-type: none"> • 90% LOD: 180 micrograms is the lowest amount that the device can detect. • 50% LOD: 100 micrograms is the lowest amount that the device can detect. • 10% LOD: 20 micrograms is the lowest amount that the device can detect. To enable a semi-quantitative device that indicates whether a misoprostol tablet is of good quality (using the cut-off value of 180 micrograms), the LOD should be set at 90%.
22B: Limit of Quantification (LOQ)	Not applicable for semi-quantitative devices	Not applicable for semi-quantitative devices	*Only applicable for quantitative tests


Variable	Minimum	Optimum	Annotations
23B: Sensitivity and specificity for detecting the target substance at the limit of detection	Sensitivity 80% Specificity 80%	Sensitivity 95% Specificity 95%	High sensitivity is crucial to minimize false negatives, ensuring that substandard or falsified drugs are detected. High specificity is essential to minimize false positives, ensuring that genuine products are not incorrectly identified as substandard.
24B: Calibration	Calibration of the device is needed, can be performed by the user Calibration time is maximum 5 minutes	Calibration of the device is not needed	
25B: Robustness	Method should be robust to different analytical procedure parameters, ambient temperature, humidity and duration of the procedure.	Same as minimum	
26B: Reproducibility	Method should give reproducible results when device is used by a different person, on a different day and in a different place after transport.	Same as minimum	
27B: Maintenance Requirements	LMIC-based third party vendors can provide maintenance. Spare parts are available within 4 weeks after placement of the order Recalibration or servicing are needed every two years	Maintenance can be performed by user itself Spare parts are available within 2 weeks after placement of the order Recalibration or servicing are needed every three years	Devices requiring minimal maintenance but with accessible support pathways are essential for long-term use in resource-limited settings. Reliance on international service may not be sustainable or responsive.

Variable	Minimum	Optimum	Annotations
28B: Internal Control/ Reference	Optional internal control/ reference items (e.g. reference library, values, standards) are either included or inexpensive and readily available in LMIC.	No internal control needed	Built-in controls improve reliability, ensure test validity, and enhance user trust— particularly critical when devices are used in decentralized settings.
29B: Companion equipment	Small, portable, reusable handheld device.	No companion equipment required.	Eliminating reliance on external equipment simplifies implementation and improves usability in clinics with limited infrastructure and in home settings.
Cost & Financial Sustainability			
31: Device cost	Single use device: Under \$5 Multi use device: Under \$1000	Single use device: Under \$2 Multi use device: Under \$100	Upfront cost of the device/technology. Should be priced appropriately for intended markets, or supplier should be willing to negotiate on pricing.
32: Cost per test	If consumables and/or single use items (for example sample holder, cartridge) are required, one test should cost less than \$2.	No consumables and/or single use items are required, there are no additional costs per test	*Excluding costs for personnel
33: Cost of software	Software is offered free of charge and made available open-source.	Same as minimum	*Only applicable if software is required

Variable	Minimum	Optimum	Annotations
34: Accessibility	Company provides shipment to every country with appropriate documentation for import of device and consumables.	Company provides shipment free of charge to every country with appropriate documentation for import of device and consumables.	
35: Availability	Robust supply chain management for devices and consumables with minimized risk of stockouts in case of increased demand. Components to operate the device are readily available from several LMIC sources.	Same as minimum	
36: Volume estimates	Volumes aligned with selective deployment of testing tools at key service delivery points in high-need areas. This scenario assumes adoption by district hospitals, referral centers, and selected community health centers that routinely provide medical abortion. The goal is to ensure targeted quality assurance where misoprostol is most frequently used, without requiring broad coverage across all facilities or pharmacies.	Volumes aligned with widespread use of testing tools across frontline health service networks and retail supply points. Used for several indications including medical abortion, PPH, and gastric ulcer management. This scenario assumes routine integration in district and sub-district hospitals, community health centers, midwife-led clinics, and pharmacy chains (public and private)	For higher-end analytical platforms (e.g., near-infrared spectroscopy), broader interest and purchasing willingness are expected, as these devices could later be extended to test other APIs beyond misoprostol—such as mifepristone, oxytocin, or other essential medicines—once methods are developed. While this may not directly increase initial testing volumes, it enhances value justification and procurement rationale for investing in a versatile, multi-analyte testing system.
37: Stability / Shelf Life / Life cycle	Single use device: Shelf life ≥ 1 years at 30–45°C and 75–90% RH.	Single use device: Shelf life ≥ 2 years at 30–45°C and 75–90% RH.	For supply chain reliability, products should remain functional across broad temperature/humidity conditions typical of

Variable	Minimum	Optimum	Annotations
	<p>Multi use device: The lifespan of the device is approximately 5 years, depending on usage intensity, environmental conditions, and maintenance</p>	<p>Multi use device: The lifespan of the device is approximately 10 years, depending on usage intensity, environmental conditions, and maintenance</p>	<p>LMICs. Shelf stability without cold chain simplifies logistics and minimizes stockouts.</p> <p>Computers and tablets, potential companion devices, usually start showing issues after 4-5 years, QC instruments have a longer life cycle but need to be replaced after a certain period of use</p>

163



3 References

1. Torloni, M. R., Bonet, M., Betrán, A. P., Ribeiro-Do-Valle, C. C. & Widmer, M. Quality of medicines for life-threatening pregnancy complications in low- And middleincome countries: A systematic review. *PLoS ONE* vol. 15 (2020).
2. Bizimana, T., Hagen, N., Gnegel, G., Kayumba, P. C. & Heide, L. Quality of oxytocin and misoprostol in health facilities of Rwanda. *PLoS One* (2021) doi:10.1371/journal.pone.0245054.
3. Bower, J. *et al.* Quality testing of mifepristone and misoprostol in 11 countries. *Int. J. Gynecol. Obstet.* (2024) doi:10.1002/ijgo.15148.
4. Hagen, N., Khuluza, F. & Heide, L. Quality, availability and storage conditions of oxytocin and misoprostol in Malawi. *BMC Pregnancy Childbirth* (2020) doi:10.1186/s12884-020-2810-9.
5. Caillet, C. *et al.* A comparative field evaluation of six medicine quality screening devices in Laos. *PLoS Negl. Trop. Dis.* (2021) doi:10.1371/journal.pntd.0009674.
6. Guo, S. *et al.* Pilot testing of dipsticks as point-of-care assays for rapid diagnosis of poor-quality artemisinin drugs in endemic settings. *Trop. Med. Health* (2016) doi:10.1186/s41182-016-0015-8.
7. Ferreyra, C. *et al.* Developing target product profiles for *Neisseria gonorrhoeae* diagnostics in the context of antimicrobial resistance: An expert consensus. *PLoS One* (2020) doi:10.1371/journal.pone.0237424.
8. Lewin, S. R. *et al.* Multi-stakeholder consensus on a target product profile for an HIV cure. *The Lancet HIV* (2021) doi:10.1016/S2352-3018(20)30234-4.
9. Pellé, K. G. *et al.* Electronic clinical decision support algorithms incorporating point-of-care diagnostic tests in low-resource settings: A target product profile. *BMJ Glob. Heal.* (2020) doi:10.1136/bmjgh-2019-002067.
10. World Health Organization [WHO]. *Target product profiles for needed antibacterial agents: enteric fever, gonorrhea, neonatal sepsis, urinary tract infections and meeting report.* <https://iris.who.int/bitstream/handle/10665/332124/9789240003897-eng.pdf?sequence=1> (2020).
11. World Health Organization. Target product profile directory. <https://www.who.int/tools/target-product-profile-database>.
12. The Food and Drug Administration - Thailand. *The Guideline for Preparation of Documents and Manual for Re- Classification of Human Modern Drugs (draft).* (2017).
13. Tyndall, A., Du, W. & Breder, C. D. Regulatory watch: The target product profile as a tool for regulatory communication: Advantageous but underused. *Nature Reviews Drug Discovery* (2017) doi:10.1038/nrd.2016.264.
14. Murtagh, M., Blondeel, K., Peeling, R. W., Kiarie, J. & Toskin, I. The relevance of target product profiles for manufacturers, experiences from the World Health Organization initiative for point-of-care testing for sexually transmitted infections.

- Arch. Public Heal.* (2021) doi:10.1186/s13690-021-00708-y.
15. United States Pharmacopoeial Convention. USP Technology Review: Global Pharma Health Fund–Minilab™. (2020).
 16. Health and Safety Executive. *Manual handling at work: A brief guide.* <https://www.hse.gov.uk/pubns/indg143.PDF>.
 17. ThermoFisher Scientific. TruScan™ RM Handheld Raman Analyzer. <https://www.thermofisher.com/order/catalog/product/TRUSCANRM>.
 18. Russek-Cohen, E., Feldblyum, T., Whitaker, K. B. & Hojvat, S. FDA perspectives on diagnostic device clinical studies for respiratory infections. *Clin. Infect. Dis.* (2011) doi:10.1093/cid/cir056.

DRAFT