

WHO International Scheme to Evaluate Household Water Treatment
Technologies

**General Testing Protocol #10:
Solar (UV and heat) Batch Systems (with or without pre-filtration
and/or disinfection addition)**

Geneva, Switzerland

September 2014

Table of Contents

1.	Product Information Requirements	1
2.	Purpose.....	1
3.	Method	2
3.1.	Replicate samples	2
3.2.	Test waters.....	2
3.3.	Microbiological Organisms and Challenge Concentrations.....	4
3.4.	Other Test Details.....	8
3.5.	Acceptable reduction deviation	12
3.6.	Records.....	12
3.7.	Completeness.....	12
4.	Procedure.....	13
5.	Defiinitions.....	16
6.	References	19

List of Tables

Table 1	General Test Water Characteristics	3
Table 2	Challenge Test Water Characteristics	4
Table 3	Microbiological Organisms and Reduction Requirements	5
Table 4	Sampling Schedule	13

1. PRODUCT INFORMATION REQUIREMENTS

Manufacturer is to provide detailed product information as required in the Expression of Interest (EOI) which is located on the WHO website at: http://www.who.int/household_water/scheme/en/. This information is to include the basic information necessary to identify the product and conduct the testing according to the manufacturer's use instruction, which may include, but not limited to:

Flowing system (gravity fed batch stand alone):

Flow rate

Volumetric capacity

Power requirements

Operating pressure

Maximum operating pressure

Operation instructions – to include: assembly, conditioning, and use instructions, daily operation and maintenance, replacement components, cleaning, backwashing and short term storage instructions (if any).

Manufacturer capacity and supporting information upon which capacity is based.

Designated Test Laboratory

The designated testing laboratory shall be identified.

2. PURPOSE

The household water treatment (HWT) product shall be evaluated for microbiological performance based on recommendations and testing principles set forth in the World Health Organization's *Evaluating Household Water Treatment Options: Health-based targets and microbiological performance specifications* (2011). Testing conducted by one of the WHO designated testing laboratories shall also be done in line with the terms and conditions outlined in the WHO International Scheme to Evaluate Household Water Treatment Technologies Procedure ("Procedure"). The Procedure can be found on the aforementioned WHO website.

The evaluation approach is based on the expected environmental conditions of the SODIS region are validated only for a specified region, this being between 35°N and 35°S. The most favourable regions for SODIS are located between latitude 15°N and 35°N and 15°S and 35°S. These semi-arid regions are characterized by the highest amount of solar radiation. Over 90% of the sunlight reaches the earth at 15° and 35° as direct radiation when there is limited cloud cover and rainfall with less than 250mm rain and usually more than 3000 hours of sunshine annually. The second most favourable region lies between the equator and latitude 15°N and 15°S. Due to high humidity and frequent cloud cover, however, the amount of scattered radiation in this region is high (about 2500 hours of sunshine annually). It is important to note that the majority of developing regions are located between latitudes 35°N and 35°S and therefore may rely on solar radiation as an energy source for solar disinfection of drinking water.

3. METHOD

3.1. Replicate samples

Three (3) production units shall be selected and run as triplicates (3) in two (2) test waters. When the batch size of the system is less than the volume required for sample collection, six (6) production units shall be required.

Systems are required to have a mechanism that alerts the user that treatment is complete. The indicator may operate on any single component of solar technology, such as heat or UV or both.

For testing purposes, containers provided with the system shall be used. If containers are not provided with the system as sold, the containers used for testing shall be based on the direction provided with the product use instructions. In the absence of direction on containers, PET containers no larger than 3L with a width (depth of water perpendicular to the sun) that does not exceed 10cm shall be used. Manufacturer are to include a statement of the containers to not exceed 10cm or data evidence that supports a greater width is acceptable and clarification of the conditions, such as a longer exposure time. The manufacturer shall provide the (new) containers for testing.

3.2. Test waters

Test water shall be prepared daily. An important aspect is that testing will be simulated to model actual field and use conditions. Two (2) a general test water (GTW) representing high quality groundwater or rainwater and a challenge test water (CTW) with more aggressive water specifications to representing surface-water. The GTW is not technology specific, and is the same for all products. The CTW, however, is based on the product's technology. Tables 1 and 2 provide the required test water characteristics and adjustment materials for all technologies. Following test water preparation, total residual chlorine, pH, turbidity, temperature, total dissolved solids (TDS), and alkalinity shall be measured and reported on the test water tank. TOC is verified during test water prep as the weight of the adjustment material addition. Sufficient volume shall be collected to allow for a retain volume for back-up analysis, if needed. The following methods, or equivalent, shall be used:

- Chlorine (total): SM 4500-Cl G or UNE-EN ISO 7393-1
- pH: SM 4500 H+ B
- Turbidity: EPA 180.1
- Temperature: SM 2550
- TDS: SM 2540C
- Alkalinity: SM 2320-B
- TOC: humic or tannic acid addition to the test water volume is to be weighted out based on the carbon content of the humic or tannic acid and is calculated to be within the test water specification range. As an alternate, SM 5310C, in water (GTW, lower TOC); SM 5310B, in water (CTW, higher TOC) may be used.

3.2.1. General test water

The general test water represents non-stressed phase of testing. Reverse osmosis treated water shall be used as the base water and adjusted to meet the following characteristics:

Table 1: General Test Water Characteristics

Constituent	Specification	Adjustment Materials (CAS# ³)
Chlorine ¹ (mg/L)	< 0.05	None
pH	7.0 ± 0.5	Inorganic acid or base: Hydrochloric acid (7647-01-0) Sodium hydroxide (1310-73-2)
TOC (mg/L)	1.05 ± 0.95 mg/L	Tannic acid (1401-55-4, Supplier: Alfa Aesar)
Turbidity (NTU)	< 1 NTU	No adjustment
Temp (°C)	20 ± 3°C	Not applicable
TDS (mg/L)	275 ± 225 mg/L	Sea Salts, Sigma Chemical Company (7732-18-5)
Alkalinity ² (mg/L as CaCO ₃)	100 ± 20 mg/L	Sodium bicarbonate (144-55-8)

¹ All chlorine shall be removed to below detection limits without the aid of added chemical(s) and is commonly accomplished by using activated carbon. Chlorine shall be measured prior to addition of test water adjustment materials. Chloride levels in Challenge Water may cause interference with analytical technique; measurements shall be made prior to addition of sodium chloride.

² Intended to buffer pH. Analyzed values may deviate from this range.

³ Chemical Abstract Service registration number. Refer to the definition section of this document for additional information.

3.2.2. Challenge Test Water

The CTW is intended for the stressed challenge phase of testing. Reverse osmosis treated water shall be used as the base water and adjusted to meet the following characteristics:

Table 2: Challenge Test Water Characteristics

Constituent	Specification	Adjustment Materials (CAS # ²)
Chlorine ¹ (mg/L)	< 0.05	None
pH	7.0 ± 0.5	Inorganic acid or base: Hydrochloric acid (7647-01-0) Sodium hydroxide (1310-73-2)
TOC (mg/L) ³	15 ± 5 mg/L	Humic acid (6813-04-4, Supplier: Alfa Aesar)
Turbidity (NTU) ³	40 ± 10 NTU	ISO spec. 12103-A2 fine test dust
Temp (°C)	4 ± 1°C ⁴	Not applicable
TDS (mg/L)	1500 ± 150 mg/L	Sea Salts, Sigma Chemical Company (7732-18-5)
Alkalinity ⁴ (mg/L as CaCO ₃)	100 ± 20 mg/L	Sodium bicarbonate (144-55-8)

¹ All chlorine shall be removed to below detection limits without the aid of added chemical(s) and measured prior to addition of test water adjustment materials) and is commonly accomplished by using activated carbon. Chlorine shall be measured prior to addition of test water adjustment materials. Chloride levels in Challenge Water may cause interference with analytical technique; measurements shall be made prior to addition of sea salts.²Chemical Abstract Service registration number. Refer to the definitions of this document for additional information.

³TOC and Turbidity added only at microbiological challenge points, except during a ‘clogging point’ during which all test water may have elevated TOC and turbidity, depending on the product specific test plan.

⁴ If product literature indicates product use limited to warm climates, the cold specification shall not be required and the GTW temperature or product appropriate temperature shall be used

⁵Intended to buffer pH. Analyzed values may deviate from this range.

3.3. Microbiological Organisms and Challenge Concentrations

Table 3 shows the organisms and American Type Culture Collection numbers (ATCC) used in evaluating performance for all technologies. The target pre-treatment concentrations of the organisms for all technologies shall be sufficient to demonstrate: *highly protective, protective, or limited protection*.

Table 3: Microbiological Organisms and Reduction Requirements

Organism	Pre-treatment Challenge ¹	Minimum Required Reduction (log)	
		Highly Protective	Protective or Limited Protection
Bacteria: <i>E. coli</i> (ATCC 11229)	$\geq 10^5/100$ mL	≥ 4	≥ 2
Virus ^{2,3} : MS-2 coliphage (ATCC 15597-B1, with host organisms: <i>E. coli</i> (ATCC 15597) or <i>Salmonella typhimurium</i> (WG49 NCTC 12484) and phiX-174 coliphage (ATCC 13706-B1) with host organisms: <i>E. coli</i> (ATCC 13706 or ATCC 700078)	$\geq 10^8/L$	≥ 5	≥ 3
Cyst: <i>Cryptosporidium parvum</i> infectious oocysts	$\geq 5 \times 10^5/L$	≥ 4	≥ 2

¹ The pre-treatment challenges may constitute greater concentrations than would be anticipated in source waters, but these are necessary to properly test, analyze, and quantitatively determine the indicated log reductions. The pre-treatment challenge must not be less than that required to demonstrate the geometric mean and standard deviation minimum required reduction described below.

² Virus performance claim will be based on the poorest log reduction of the two phages.

³ Host selection is dependent on method. Refer to Section 3.3.1.4 Organism Methods.

3.3.3. Selection of Microorganisms

It is not practical, and there are insufficient data, to set performance targets for all potentially waterborne pathogens. Therefore, the most sensible approach is to identify reference pathogens that represent groups of pathogens. The Scheme reference target organisms were chosen to represent classes of pathogens in water (bacteria, virus and protozoa) with respect to occurrence, concentration and health impact.

For actual testing of performance, selection of microorganisms that represent the three classes of pathogens is necessary. Ideally, surrogates would be chosen for all classes as they are easier and cheaper to use, two important considerations for making the protocol accessible to range of laboratories. However, at this time, there is insufficient data to support selecting surrogates for all classes of pathogens. Thus, the microorganisms selected for inclusion for the HWT Scheme are well documented as laboratory test organisms; they have varying degrees of susceptibility to commonly used drinking water disinfectants; and represent an array of particle sizes/surface properties that should provide useful information with respect to HWTs that rely on mechanical size exclusion for the reduction of microbes.

3.3.3.1. Enteric Bacteria

Enteric bacteria are generally the group of pathogens most sensitive to inactivation by disinfection. The bacteria species *Escherichia coli* (*E. coli*) shall be used to represent the challenge of bacterial contaminants. *E. coli*, as well as members of the *Enterobacteriaceae* family, has a history of use in disinfection studies and protocols. *E. coli* is typical of the total coliform bacteria group frequently found in untreated surface waters and has added

health significance as its presence is very indicative of fecal contamination. Some strains of *E. coli* produce toxin(s) that can lead to severe gastrointestinal illness. According to a recent global study in over 20,000 children in seven developing countries, *E. coli* was among the top three pathogens associated with moderate to severe diarrhoea (Kotloff, et al., 2013).

3.3.3.2. Enteric Virus

Human enteric viruses are the smallest pathogens, making them more difficult to remove by physical processes, such as filtration. Specific viruses may be less sensitive to disinfection than enteric bacteria and some protozoan parasites. Using human or animal viruses in laboratory testing is complicated, expensive and given the availability of comparable surrogates, this latter option was chosen for the Scheme. Two different surrogate bacteriophages, MS-2 and phiX-174, shall be used to evaluate the performance of HWT products for performance. In choosing surrogates, consideration included the wide variety of different viruses' resistance to potential treatment processes that enteric viruses vary greatly in terms of size, isoelectric points, type of nucleic acid, presence of lipids, and the structure of the proteins in the capsid. Additionally, some treatment systems have more than one mechanism that would remove/inactivate viruses. For example, a filtration system (activated carbon) may be combined with a UV light system. Some viruses may be more easily removed by adsorption to the activated carbon than others, and others may be more resistant to the UV light. For these reasons and due to not using an actual pathogen, the testing of two bacteriophages, with varying characteristics and responses to treatment processes, shall be used in the assessment of the performance of HWT products.

MS-2 and phiX-174 are extensively used bacteriophages as models for human enteric virus removal by water treatment processes. A great deal is known about the resistance of these bacteriophages to disinfectants. They are easy to grow to large number. Both are similar size and lack a lipid coat like many of the human enteric viruses.

MS-2, 24nm in diameter, is a singled stranded RNA virus, with a low isoelectric point (3.9). A low isoelectric point makes it less sticky (adsorbs to a lesser degree) than poliovirus and has been used as a conservative model for removal by adsorption processes. MS-2 is very resistant to inactivation by low-pressure UV light and has been used as a model virus to measure UV light dose in UV light reactors (collimated beam). It is one of the more hydrophobic non-lipid containing viruses.

PhiX-174, 25-27nm in diameter, is a single stranded DNA with an isoelectric point of 6.6. It is less hydrophobic than MS-2. Research suggests that it is more resistant to halogen disinfectants like iodine and chlorine dioxide than MS-2.

3.3.3.3. Parasitic Protozoa

The oo/cysts of parasitic protozoa are the group of pathogens least sensitive to inactivation by chemical disinfection, but relatively sensitive to UV light irradiation, as seen with oocysts of *Cryptosporidium*, which are highly resistant to oxidizing disinfectants such as chlorine. Protozoan oo/cysts are of a moderate size (>2um) and are more readily removed by physical processes compared to viruses and bacteria. Causing the disease Cryptosporidiosis, a severe gastrointestinal illness, *Cryptosporidium hominis* and *C. parvum* are pathogens of concern worldwide and key waterborne reference pathogens cited in the GDWQ (WHO, 2011). According to the same recent, aforementioned study, *Cryptosporidium* is one of the top three pathogens responsible for diarrhoea in young children in developing countries (Kotloff, et al., 2013). In the environment, the organism exists in a protective cyst stage called an oocyst. *Cryptosporidium* oocysts are typically 3-5 microns in diameter, making it a suitable representative to challenge filtration technologies. *Cryptosporidium parvum* infectious

oocysts shall represent the challenge to evaluation protozoa reduction and/or inactivation performance.

For all testing, a total of 1.4L of product water shall be collected and sub-sampled based on analysis sample volume requirements for each microbiological test organism. The 1.4L is sufficient sample size to allow for organism analysis and a retain volume.

3.3.3.4. Organism methods

Production and assay procedures for the microbial challenges and equivalent methods shall include, but not be limited to:

- ***E. coli*** (ATCC 11229) shall be prepared using the method specified in Asburg, E.D. Methods of Testing Sanitizers and Bacteriostatic Substances; in *Disinfection, Sterilization, and Preservation* (Seymour S. Block, ed.) (1983). The samples shall be assayed in triplicate with m-Endo medium using Method 9222B in Standard Methods for the Examination of Water and Wastewater (APHA, 2012). The geometric mean and standard deviation of the triplicate assay shall be reported for each water type and across all water types examined.
 - Collected samples shall be stored at a temperature between 1 – 8°C and processed within 24 hours.
 - Required sample volume to allow for processing in triplicate and a retain volume: 660ml
- **Coliphage MS-2** (ATCC 15597-B1) shall be prepared and assayed using:
 - The method in Annex A, Section A.8.2.2 of NSF/ANSI 55: Ultraviolet Microbiological Water Treatment Systems (2012); *E. coli* host ATCC 15597; or
 - NEN-EN-ISO 10705-1 (Detection and enumeration of bacteriophages Part 1: Enumeration of F-specific RNA bacteriophage).
 - *Salmonella typhimurium* (WG49) host NCTC 12484 or *E.coli* host ATCC 15597. Analyses shall be conducted in triplicate; the geometric mean and standard deviation for each water type and across all water types examined shall be reported.
 - Samples shall be stored at a temperature between 1 – 8°C and processed within 24 hours of collection.
 - Required sample volume to allow for processing in triplicate and a retain volume: 12ml
- **Coliphage phiX-174** (ATCC 13706-B1) shall be prepared and assayed using:
 - The method in Annex A, Section A.8.2.2 of NSF/ANSI 55: Ultraviolet Microbiological Water Treatment Systems (2012); *E. coli* (host) ATCC 700078; or
 - NEN-EN-ISO 10705-1 (Detection and enumeration of bacteriophages Part2: Enumeration of somatic coliphages)
 - *E. coli* host ATCC 700078 or ATCC 13706
 - Analyses shall be conducted in triplicate; the geometric mean and standard deviation for each water type and across all water types examined shall be reported.
 - Samples shall be stored at a temperature between 1 – 8°C and processed within 24 hours of collection.

- Required sample volume to allow for processing in triplicate and a retain volume: 12ml
- *Cryptosporidium parvum* infectious oocysts shall be assayed using an infectivity method which shall be based on a “Most-Probable-Number Assay (MPN) for Enumeration of Infectious *Cryptosporidium parvum* Oocysts”, including the standard deviation, as per Slifko *et al.* (1999) for each water type and across all water types examined.
 - Samples shall be stored at a temperature between 1 – 8°C and concentrated by centrifugation within 24 hours of collection. R
 - Required sample volume to allow for processing in triplicate and a retain volume: 600ml

3.4. Other Test Details

3.4.4. Solar simulation

Emissions from the sun include visible light, heat (infra-red) and UV radiation. The Atlas Suntest XXL with the ‘coated quartz/daylight’ flat optical filter installed shall be used to provide simulated outdoor solar radiation. This configuration shall deliver an irradiant range of 300nm – 800nm, which is equivalent to 550W/m², and provides infra-red, visible and UV-A.

3.4.5. Untreated control

The microbiologically spiked test water to be used as the pre-treatment challenge concentration, shall also serve as the untreated control. See Table 3 for concentrations.

3.4.6. Blank sample

Prior to test initiation, using the GTW, the product shall be tested for the presence of the test organisms without microbiological addition to confirm that the product arrived to the laboratory free of test organisms. Systems shall flow sufficient volume of GTW, with no microbiological addition, through the system challenge to allow for the collection of the necessary volume for analysis for the organisms of Table 3.

3.4.7. Quality assurance/quality control (QA/QC)

The testing laboratory will adhere to the requirements of their QA/QC procedures and ISO 17025 requirements and must be able to provide documentation of adherence, which are to include but not be limited to quality checks on organism stocks, calibration of instruments, testing environmental controls, etc.

3.4.8. Product residual or wetted contact material of concern

For products that employ a disinfectant, bacteriostatic agent, or have a wetted contact material which may have a contaminant leach concern, product residual sample shall be collected with the microbiology samples at each microbiological challenge point. The active agent residual shall not constitute a threat to health. The WHO Guidelines for Drinking-water Quality (2011) shall be used to determine acceptable levels in the product water.

3.4.9. Neutralization

For products that employ a disinfectant, verification of the efficacy of neutralization of the product residual shall be verified for both test waters (GTW and CTW). The Untreated Control shall address potential issues of toxicity of the neutralizer. The methods are described in ASTM E1054-08 (2013).

3.4.10. Microbiological sample points

Sampling for microbial pathogens shall be conducted according to the schedule shown in Table 4. Seeding shall be required for microbiological challenge points for systems which have pre-filtration of other components which hold a volume between batches.

- Seeding shall be used to purge the system of the uncontaminated water with a sufficient flow of contaminated test water (seeding). The systems shall be exposed to a minimum of 10 units void volumes or 1L, whichever is greater, of microbiologically challenged water per Table 3 immediately prior to sample collection and continued through sample collection.
- A 'sampling event' includes seeding and sample collection.
- For batch systems, a full batch shall be used for seeding and a full batch shall be collected and sub-sampled into prepped bottles for microbiological analysis.
- If a product requires seeding, based on expected length of time to treat a batch, it may be necessary to seed the system on one day and treat a sample for collection the following day. The laboratory shall have the option to conduct seeding and batch treatment for sample collection consecutively, allowing 6.5 hours per batch for treatment according to Section 3.4.10.
- If the system does not have any pre-filtration or other components that will hold a volume between batches, seeding would not be required. A full batch shall be collected and sub-sampled into prepped bottles for microbiological analysis.

Sample Collection for GTW:

- 1 blank sample analyzed for the organism of Table 3 and any bacteriostatic or disinfection chemical and any wetted material contaminant of concern, if system includes such.
- 3 pre-treatment sample analyzed for organism of Table 3
- 1 sample for Test Water Characteristics
- 9 post-treatment samples analyzed for organism of Table 3 and the bacteriostatic or disinfection chemical and any wetted material contaminant of concern, if system includes such.

Sample Collection for CTW:

- 1 blank sample analyzed for the organism of Table 3
- 3 pre-treatment sample analyzed for organism of Table 3
- 1 sample for Test Water Characteristics
- 9 post-treatment samples analyzed for organism of Table 3, any bacteriostatic or disinfection chemical and any wetted material contaminant of concern, if system includes such.

Flowing systems with chemical disinfection shall require neutralization and residual disinfection concentration analysis as discussed under Section 3.4.6.

3.4.11. Conditioning

Conditioning, if required, shall be according to the Operation Manual and shall use GTW without microbiological addition. Conditioning volume shall not accumulate against the test capacity volume.

3.4.12. End of life

The manufacturer must provide an explicit indication or assurance of the unit's effective use lifetime to warn the consumer of potential diminished treatment capacity by one of the following:

- Having the unit terminate discharge of treated water
- Sounding an alarm
- Providing single explicit instructions for servicing or replacing units within the recommended use life (measurable in terms of volume throughput, specific timeframe or other appropriate method).
- Or a statement of the understanding that the system will clog.

3.4.13. Daily test capacity

The manufacturer provided use instructions/operator's manual may supply the daily capacity of the system and the system shall be run accordingly, but not to exceed 6.5 total hours of system operation in a single test day. The daily test capacity will be based on product use, time for treatment and laboratory efficiency.

3.4.14. Leakage test

Flowing systems shall not leak during test operation. Any leaking during test operation shall be recorded in the laboratory bench sheets.

3.4.15. Device cleaning

Approaches to restore or maintain flow identified in the user instruction/operator manual shall be permitted during testing.

- The test plan shall detail the trigger to initiate a cleaning or backwash, such as 75% flow reduction, or may be indicated to be performed routinely after an identified accumulated volume, such as every 10L or prior to storage.
- Device cleaning shall not occur during seeding and sample collection at any sample point, even if the trigger point is achieved during the seeding or collection of sample.
- Changes in flow rate and cleaning performed during the testing shall be recorded in the laboratory bench sheets.
- Unless the product literature includes information that would identify a different product specific trigger, 75% reduction shall be identified in the test plan.

3.4.16. Component replacement

For systems, a component that would not be considered a primary component in providing the microbiological reduction performance may be replaced as needed during the test. An example is a pre-filter for turbidity removal. However, a component which provides microbiological performance shall not be replaced during the testing. The general test plan for the specific technology provides direction on component replacement during testing.

3.4.17. End of test

The general test plan for the product type shall provide clear direction on ‘end of test’. For flowing system devices, there shall be three (3) acceptable outcomes for the end of the test:

1. Completion of the Sampling Schedule of Table 4.
 - a. The test is complete after the collection of the 6 test days (test capacity) according to Table 4.
 - b. The test capacity (volume) shall be the accumulated volume during the 6 day test.
2. The failure of more than one replicate to indicate treatment complete during the daily test capacity as defined in Section 3.4.10 for more than a single test day.
 - a. If a replicate’s indicator does not trigger for treatment complete, there shall be no sample collected from the replicate.
3. Clogging, as a reduction of greater than 75% flow compared to initial flow that cannot be restored with test permitted component replacement or cleaning procedures.
 - a. The test capacity, as volume, shall be the accumulated volume during the time on test.

3.4.18. Log reduction calculation

Testing shall be conducted simultaneously on the technology dictated number of replicates. At each microbiological sampling point, pre-treatment and post-treatment water samples shall be collected and each analyzed in triplicate.

When reporting the geometric means of the triplicate counts, if all three counts are non-detect for the organism, the geometric mean should be reported to indicate “Less than” (<). In the event one or more PFU, CFU, or oocysts are found in one or two of the triplicate counts, the “less than” counts are to be treated as being at the detection limit for the purpose of calculating the geometric mean and standard deviation.

Log reductions for the purpose of compliance with this test plan shall be calculated at each sample point as follows:

The geometric mean (GM) of each triplicate analysis (X) shall be calculated for each pre-treatment sample and replicate post-treatment as:

$$[1] \quad GM = (X_1 * X_2 * \dots * X_n)^{(1/n)}$$

The geometric mean is defined as the n th root (where n is the count of numbers) of the product of the numbers. Such as, the geometric mean of the three numbers is the cube root of their product.

The geometric mean applies only to positive numbers. It is also often used for a set of numbers whose values are meant to be multiplied together or are exponential in nature, such as data that will be reported for the microbiological concentration in the pre-treatment and post-treatment waters of the testing.

Example of calculating the geometric mean:

Use triplicate post-treatment analyses results of: 1.00E+02, 7.70E+01, and 9.30E+01.

Since there are 3 numbers, the n -th root is the 3rd root. The geometric mean would be:

$$(1.00E+02 * 7.70E+01 * 9.30E+01)^{1/3} = 8.95E+01$$

The log reduction for each replicate at each sample point shall be calculated using the results from [1], shown below as the negative \log_{10} of the GM of each replicate post-treatment, GM_{eff} , divided by the GM of the pre-treatment, GM_{inf} .

$$[2] \quad \text{Log Reduction} = -\log_{10}(GM_{\text{eff}}/GM_{\text{inf}})$$

Example of calculating the log reduction:

Using an example pre-treatment geometric mean of 2.07×10^8 units (such as CFU/100mL) and using the above examples reported geometric mean of the triplicate analysis of 8.95×10^1 units (CFU/100mL), the log reduction would be:

$$\begin{aligned} \log_{10} 2.07 \times 10^8 - \log_{10} 8.95 \times 10^1 &= \text{log reduction} \\ 8.31 - 1.95 &= 6.36 \text{ log reduction} \end{aligned}$$

For reporting purposes, two (2) significant figures shall be reported. For the above example, 6.4 would be reported. For evaluation of log reduction against the pass/fail criteria, ASTM Standard E29 Absolute method shall be used, which does not allow for rounding.

3.5. Acceptable reduction deviation

Three (3) production products for three (3) lots of the product must continuously meet or exceed the reduction requirements shown in Table 3, except for the following acceptable allowance. Up to 10% of pre-treatment/post-treatment sample pairs may vary from the reductions required in Table 3 but not achieve less than:

Viruses:	1 log variance
Bacteria:	1 log variance
Oocysts:	1 log variance

Each phage is treated separately for evaluating acceptable allowance, however the overall claim for virus shall be based on the phage for which the product performed the poorest. Additionally, the geometric mean of all microbiological reductions must meet or exceed the requirements of Table 3. Compliance with the requirements shall be based on the reduction percentage calculation.

3.6. Records

All pertinent procedures and data shall be recorded and provided in a final report. The general test plan for the product type provides a list of the data that is to be reported.

3.7. Completeness

Completeness is a measure of the number of valid samples and measurements that are obtained during a test period. Completeness will be measured by tracking the number of valid data results against the specified requirements in the test plan.

Completeness will be calculated by the following equation:

$$\text{Percent Completeness} = (V / T) \times 100\%$$

Where:

V = number of measurements that are valid

T = total number of measurements planned in the test

The specification for this data quality objective will be to achieve minimum 90% completeness for microbiological and disinfectant residual samples scheduled in the test plan or one (1) incomplete measurement (if less than 10 are taken).

4. PROCEDURE

Test waters shall be prepared daily and verified in accordance with Tables 1 and 2, except the conditioning water. Daily test water characteristics shall be sampled, analyzed and results provided in the final report. All sample volumes collection, both microbiological and chemical shall be collected such that sufficient sample volume remains after analysis to allow for retain sample. The remaining volume of sample shall be retained for confirmation or retesting purposes, when necessary. Flows shall be measured and recording at each test day start, at sample collection (start), and at the end of each test day.

Table 4: Sampling Schedule

Test Day ¹ : Collection point	Test Water	Microbiological Sampling	
Day 1: At indication of “treatment complete”	GTW	Pre-treatment	Post-treatment
Day 2: At indication of “treatment complete”	GTW	Pre-treatment	Post-treatment
Day 3: At indication of “treatment complete”	GTW	Pre-treatment	Post-treatment
Day 4: At indication of “treatment complete”	CTW	Pre-treatment	Post-treatment
Day 5: At indication of “treatment complete”	CTW	Pre-treatment	Post-treatment
Day 6: At indication of “treatment complete”	CTW	Pre-treatment	Post-treatment

¹When seeding is required, each sample collection may required 2 days for seeding and sample batch processing

PROCEDURE

1. Three (3) systems or six (6) shall be selected and tested simultaneously.
2. All systems shall be conditioned according to the Operator's Manual use instructions.
3. Test waters shall be verified to be in accordance with Tables 1 and 2.
 - a. Daily test water characteristics shall be sampled, analyzed and results provided in the final report.
 - b. Sufficient volume shall be collected to allow for a retain volume.
 - i. One sample volume shall be used for analysis and reporting.
 - ii. The second volume shall be retained for confirmation or retesting purposes, when necessary.
4. The Atlas Suntest XXL with the 'coated quartz/daylight' flat optical filter installed shall be used to provide simulated outdoor solar radiation.
 - a. The light shall be applied continuously until the product's indicator displays that treatment is complete.
 - b. The daily total run time exceed 6.5 total hours of exposure in a single test day.
5. Devices shall be operated according to the Operator's Manual.
 - a. Test plan permitted system cleanings or backwashes and/or component replacements shall occur as dictated by operator's manual's requirements.
 - b. Where applicable, the laboratory technician(s) shall note treatment time, daily starting and ending flow rates, flow rates at sample collection, and any cleaning procedures, and other significant event throughout the testing.
6. Testing shall be according to the schedule of Table 4.
 - a. Devices shall be operated based on the identified daily use.
 - b. During overnight storage, the devices shall be prepped according the manufacturer's instruction for short term storage during normal operation and stored at room temperature.
 - c. Microbiological sample collection shall occur immediately following treatment complete as determined by product installed indicator, which may vary by replicate. At each microbiological challenge point, the following shall occur:
 - i. Microbiologically spiked challenge water shall be prepared to meet the concentrations of Table 3. All organisms (*E. coli*, MS2, PhiX-174 and *Cryptosporidium parvum* oocysts) are compatible for combined challenge.
When six (6) units are required, organisms pairing shall be *Cryptosporidium*/MS-2 and *E.coli*/phiX-174 (See Section 3.1).
 - ii. A complete batch of microbiologically spiked (Table 3) test water shall pass through the system as seeding, when appropriate per Section 3.4.7. If a single batch is less than 1L or does not be exposed the system to a minimum of 10 units void volumes of microbiologically challenged water, additional batches shall be passes to achieve seeding requirements.
 - iii. A full batch of microbiologically spiked test water (Table 3) shall be processed by the system, collected, and sub-sampled into appropriately prepped collection vessels as post-treatment samples for microbiological analyses.

- iv. A sample shall be collected from the pre-treatment challenge and analyzed to confirm pre-treatment concentrations.
 - 1. As the pre-treatment challenge water for all units was from a single source of microbiologically spiked prepped test water, a single pre-treatment sample shall be taken and for analyzed for the organisms of Table 3.
 - 2. The pre-treatment microbiological shall be collected immediately from the microbiologically spiked test water after all post-treatment samples have been collected.
 - v. For all sample collection, sufficient volume shall be collected to allow for a retained as a backup for confirmation or retesting purposes, when necessary.
 - vi. Microbiological pre-treatment and post-treatment concentrations shall be presented in the final report.
7. The end of test shall be one of the following:
- a. Completion of the Sampling Schedule of Table 4. The test capacity, as volume, shall be the accumulated volume during the time on test.
 - b. The failure of more than one replicate to indicate treatment complete during the daily test capacity as defined in Section 3.4.10 for more than a single test day. If a replicate's indicator does not trigger for treatment complete, there shall be no sample collected from the replicate.
 - c. Clogging, as a reduction of greater than 75% flow compared to initial flow that cannot be restored with test permitted component replacement or cleaning procedures. The test capacity, as volume, shall be the accumulated volume during the time on test.

5. DEFINITIONS

The following establishes definitions for terminology used with household water treatment as point-of-use or point-of-collection disinfectants or units and related components. This list is general for all Generic Test Plans (GTPs) established under *Evaluating Household Water Treatment Options: Health-based targets and microbiological performance specifications* (2011).

Active agent: A substance or medium added to or involved in a drinking water treatment process that requires direct or sacrificial release of the agent or its degradation product(s) to perform a specific functions.

Additive: A substance added to water, directly or indirectly, during a drinking water treatment process.

Backwash: A reversed flow of water through a media which allows the expelling of collected matter to the drain.

Back flush: The references of flow direction through a filter or ion exchange column or membrane to remove particles for cleaning purposes.

Bacteriostatic: A biological or chemical agent that stops bacteria from reproducing, while not necessarily harming them otherwise.

Batch treatment: A method in which a fixed quantity of water is processed through a treatment device in a single treatment cycle.

Capacity: The volume of water treated by a system before the system or components of the system must be cleaned, regenerated or replaced, as specified by the manufacturer.

Challenge water: The mixture of water and contaminants used to test a system for contaminant reduction claims.

Chemical Abstract Service (CAS) Registration Number (RN): Unique numerical identifiers assigned by the Chemical Abstracts Service to every chemical described in the open scientific literature (currently including those described from at least 1957 through the present) and including elements, isotopes, organic and inorganic compounds, ions, organometallics, metals, nonstructurable materials. They are referred to as CAS RNs and CAS Numbers. A CAS RN designates only one substance, has no chemical significance, and provides a link to information about a specific chemical substance. Chemical compounds can be described in many different ways such as molecular formula, structure diagram, systematic names, generic names, proprietary or trade names, or trivial names. A CAS Registry Number, however, is unique and specific to only one substance. CAS Registry Numbers allow for keeping track of substances because they are unique, can be validated quickly and reliably, and are internationally recognized. As CAS RNs are not dependent upon any system of chemical nomenclature, they can provide a reliable common link between the various nomenclature terms used to describe substances and serve as an international resource for chemical substance identifiers used by scientists, industry and regulatory bodies. The assigning agency, Chemical Abstracts Service (CAS) is a function of the American Chemical Society (ACS) and CAS information is copyrighted by the ACS. www.cas.org

Cleaning: Removal of residues and other soiling materials.

Component: A separate or distinct part of a water treatment system including, but not limited to membranes, filters, housings, tubing, storage tanks, faucets, valves, and connectors.

Oo/cyst: The environmentally resistant stage in the life cycle of certain parasitic protozoa which are identified from water samples. These include oocysts of *Cryptosporidium* and *Toxoplasma* and cysts of *Giardia* and *Entamoeba*.

Daily production rate: The volume of product water produced by the system per day under defined conditions.

Disinfection: The process that eliminates (removing, destroying, and inactivating) many or all pathogenic microorganisms with the exception of the bacterial endospore on inanimate objects and liquids.

Post-treatment: The treated water from the outlet of a unit, system, component, or process.

Filter: (verb) To pass water through a permeable medium to separate particles from the water. (noun) A device for carrying out the process of filtration consisting of the medium and suitable hardware for constraining and supporting the medium in the path of the water.

Filtration: The process by which particles are separated from water by passing water through a permeable material.

Hardness: A measurement of the concentration of divalent and trivalent cations, primarily calcium and magnesium, in drinking water. Hardness is typically expressed as grains per gallon or mg/L as calcium carbonate.

Household Water Treatment (HWT) Technology: A product that is used in households or similar settings to remove water contaminants that may pose health risks. Priority products for testing will be low-cost, appropriate for the poor and generally “free standing” products which only treat enough water to serve a limited number of individuals.

Pre-treatment challenge: The mixture of water and contaminants entering a water treatment system.

In-line device: Any device in contact with the water installed on a service line or distribution system downstream of the water main and upstream from endpoint devices.

Media: Material in a system that forms a water-permeable barrier to the passage of certain contaminants or otherwise contributes to the reduction of contaminants in water. Medium is the singular form of media.

Membrane: A semi-permeable barrier that allows the passage of water, and depending on membrane type and characteristics, may restrict the passage of microorganisms, particles, molecules, and ions.

pH: The negative log of the hydrogen ion concentration a measure of the degree of acidity or alkalinity of an aqueous solution.

Post-treatment: The treated water from the outlet of a unit, system, component, or process.

Pre-treatment challenge: The mixture of water and contaminants entering a water treatment system.

Pressure: The force applied to a unit area. Water pressure is normally measured in lb/in², kilopascals (kPa), or feet or metres of head.

Product water: Water that has been treated by the system.

Rated service cycle: The capacity or time of operation of a system or component between cleaning, replacement, or regeneration of the treatment medium (media), as specified by the manufacturer.

System: A complete water treatment device, including all components needed to connect it to a potable water supply.

Total dissolved solids (TDS): The solids remaining when a solution is filtered through a 0.45 µm glass filter and the filtrate is evaporated and dried to constant weight at 180 °C (356 °F). TDS is expressed as mg solids per litre of filtrate.

Turbidity: A condition caused by the presence of suspended matter, colloidal matter, or both, which results in the scattering and absorption of light.

Unit void volume: Total water-holding volume with the medium (media) and internal components in place.

Unit volume: Total water-holding volume without the medium (media) or internal components.

6. REFERENCES

The following references have been assembled as a single list to cover all Generic Test Plans (GTPs) established under *Evaluating Household Water Treatment Options: Health-based targets and microbiological performance specifications* (2011). As such, not all references are applicable to the GTP of this document.

American Public Health Association (APHA), American Water Works Association (AWWA), and Water Environment Federation (WEF). (2012) *Standard Methods for the Examination of Water and Wastewater*, 21st Edition.

- 4500-H⁺ pH Value, B. Electrometric Method
- 4500- Cl Chlorine (Residual), G. DPD Colorimetric Method
- 2550-Temperature
- 2540-Solids, C. Total Dissolved Solids Dried at 180°C
- 2320-Alkalinity, B. Titration Method
- 5310-Total Organic Carbon (TOC), C. Persulfate-Ultraviolet or Heated-Persulfate Oxidation Method
- 5310-Total Organic Carbon (TOC), B. High-Temperature Combustion Method
- 9222-Membrane Filter Technique for Members of the Coliform Group, B. Standard Total Coliform Membrane Filter Procedure

ASTM (2013). E1054-08. Standard Test Methods for Evaluation of Inactivators of Antimicrobial Agents. ASTM International.

ASTM (2009). C29-08. Standard Practice for Using Significant Digits in Test Data to Determine Conformance with Specifications. ASTM International.

Bingham et al. (1979). *Giardia* sp.: physical factors of excystation in vitro and excystation vs eosin exclusion as determinants of viability. Exp. Parasitol. 47:284-291.

Bingham and Meyer (1979). *Giardia* excystation can be induced in vitro in acidic solutions. Nature. 277:301-302.

Block, SS. (1983). Disinfection, Sterilization, and Preservation. Philadelphia: Lea & Febiger.

International Organization for Standardization. (1997). ISO 12103 Part 1.A2: Road vehicles - Test dust for filter evaluation, Arizona test dust, Fine Test Dust.

International Organization for Standardization. (2006). ISO 15553. Water quality – Isolation and identification of *Cryptosporidium* oocysts and *Giardia* cysts from water.

Kotloff, et al. (2013). Burden and aetiology of diarrhoeal disease in infants and young children in developing countries (the Global Enteric Multicenter Study, GEMS): a prospective, case-control study. Lancet. 382: 209–222.

NEN-EN ISO 10705-1 (1996). Detection and enumeration of bacteriophages Part2: Enumeration of somatic coliphages.

NEN-EN ISO 10705-1 (1996). Detection and enumeration of bacteriophages Part2: Enumeration of somatic coliphages.

UNE-EN ISO 7393-1 (2000). Water quality. Determination of free chlorine and total chlorine. Part 1: Titrimetric method using N,N-diethyl-1,4-phenylenediamine.

- NSF International. (2008). Protocol P231, Microbiological Water Purifiers.
- NSF International. (2012). Protocol P248, Military Operations Microbiological Water Purifiers.
- NSF International and American National Standards Institute. (2014). NSF/ANSI 53: Drinking Water Treatment Units – Health Effects. Annex A.
- NSF International and American National Standards Institute. (2014). NSF/ANSI 53: Drinking Water Treatment Units – Health Effects. Annex B.
- NSF International and American National Standards Institute. (2013). NSF/ANSI 55: Ultraviolet Microbiological Water Treatment Systems. Annex A, Section A.8.2.2.
- Rice, Hoff, Schaefer. (1982). Inactivation of *Giardia* Cysts by Chlorine. Appl Environ Microbiol. 43: 250-251.
- Rice and Schaefer. (1981). Improved in vitro excystation procedure for *Giardia lamblia* cysts. J. Clin. Microbiol. 14:709-710.
- Schaefer, et. al. (1984). Factors promoting in vitro excystation of *Giardia muris* cysts. Trans., Royal Soc. Of Trop. Med. & Hyg. 78:795-800.
- Slifko, TR, DE Huffman, JB Rose. (1999). A Most-Probable-Number Assay for enumeration of infectious *Cryptosporidium parvum* oocysts. Appl Environ Microbiol. 65(9): 3936–3941.
- U.S. Environmental Protection Agency: Office of Pesticide Programs and Office of Drinking Water. (1987). Task Force Report, Guide Standard and Protocol for Testing Microbiological Water Purifiers.
- U.S. Environmental Protection Agency. (1993) EPA Method 180.1: Determination of Turbidity by Nephelometry.
- World Health Organization. (2011). Evaluating Household Water Treatment Options: Health-based targets and microbiological performance specifications. Geneva, Switzerland. http://whqlibdoc.who.int/publications/2011/9789241548229_eng.pdf
- World Health Organization. (2011). Guidelines for Drinking-water Quality. Geneva, Switzerland. http://www.who.int/water_sanitation_health/publications/2011/dwq_guidelines/en/